

2601 Hardwood Lane
Houston, Texas 77093
1-888-668-5633



ADVANCED NUCLEAR CONSULTANTS

Philips Cardio MD Gamma Camera

The Cardio MD is a dedicated nuclear cardiac camera with dual 90° fixed detectors. Its small footprint easily fits into an 8' x 10' room and its horizontal imaging position allows for both supine and prone myocardial perfusion imaging.

Features & Benefits

- Light weight
- 4" Table cut-offs for closer imaging capabilities
- Rectangular 370mm x 214mm

ANC's Standard Equipment

- Uninterruptible Power Supply (UPS)
- Network Printer
- Cardiac IVY Gate
- Bar Phantom

ANC's Standard Services

- Equipment Applications Training
- 24 Hours on Site Response Time
- Less Than One Hour Phone Response
- Unlimited Phone Support
- Free Equipment Shipping and Freight Charges

Options to Consider

- Hot Lab and Ancillary Equipment
- Patient Comfort Kit
- Patient Table Limit of 450 pounds
- Extended Service Agreement
- Preventive Maintenance
- Archive Tracker (Data Back up)



Thank you very much for your inquiry. We appreciate the opportunity to serve you!