

2601 Hardwood Lane
Houston, Texas 77093
1-888-668-5633



ADVANCED NUCLEAR CONSULTANTS

Philips/ADAC Vertex EPIC Gamma Camera

The Philips Vertex EPIC is a digital nuclear medicine camera. The Vertex EPIC is a variable angle gamma camera utilizing large field of view dual detectors. The Vertex EPIC is designed for SPECT and whole body imaging.

Features & Benefits

- Dual ring open gantry design for stability
- Robotically controlled gantry motions
- Rectangular 200mm x 150mm FOV
- Patient table limit is 500 pounds

ANC's Standard Equipment

- Uninterruptible Power Supply (UPS)
- Network Printer
- Cardiac IVY Gate
- Bar Phantom

ANC's Standard Services

- Equipment Applications Training
- 24 Hours on Site Response Time
- Less Than One Hour Phone Response
- Unlimited Phone Support
- Free Equipment Shipping and Freight Charges

Options to Consider

- Hot Lab and Ancillary Equipment
- Patient Comfort Kit
- Extended Service Agreement
- Preventive Maintenance
- Archive Tracker (Data Back up)



Thank you very much for your inquiry. We appreciate the opportunity to serve you!