

2601 Hardwood Lane  
Houston, Texas 77093  
1-888-668-5633



# ADVANCED NUCLEAR CONSULTANTS

## Siemens Orbiter:

ANC remanufactures the Siemens Orbiter by integrating a state-of-the-art acquisition and processing system. This can be done as an upgrade to a nuclear camera currently in the field or can be purchased as a complete system directly from ANC. The remanufacturing process eliminates many of the original manufacturer's electronics, reducing potential service calls, and includes a completely new crystal assembly, a rebuilt detector and a new fast, reliable system. The unique ANC remanufacturing process results in increased performance, better diagnostics and faster imaging and data processing for all of your nuclear cardiology and general nuclear medicine procedures.

## Features & Benefits

- Improved energy, linearity and uniformity corrections
- Real-time isotope display during all acquisitions
- Intrinsic and extrinsic uniformity correction tables for any isotope-collimator pair
- Latest generation USB V2 Acquisition Interface
- Circular 387mm FOV

## ANC's Standard Equipment

- Uninterruptible Power Supply (UPS)
- Network Printer
- Cardiac IVY Gate
- Bar Phantom

## ANC's Standard Services

- Equipment Applications Training
- 24 Hours on Site Response Time
- Less Than One Hour Phone Response
- Unlimited Phone Support
- Free Equipment Shipping and Freight Charges

## Options to Consider

- Hot Lab and Ancillary Equipment
- Patient Comfort Kit
- Extended Service Agreement
- Preventive Maintenance
- Archive Tracker (Data Back up)



Thank you very much for your inquiry. We appreciate the opportunity to serve you!